



Groundwater – PLFA Sampling Protocol

SAMPLING INSTRUCTIONS

The following sampling instructions are used for collecting water or groundwater samples for PLFA analysis. The recommended sampling containers are 2-1L Poly bottles with screw caps. Amber glass bottles can be used but are not recommended due to the likelihood of breakage during shipment. Microbial Insights, Inc. can provide the proper sampling supplies upon request. Once the proper sampling bottle is obtained be sure not to contaminate the inside of the sample bottle with skin, dirt or any form of debris (this helps to ensure the accuracy of the data results). Wearing latex gloves (or similar) is recommended when sampling.

Note: It is important to collect as close to the required amounts as possible to ensure the ability to properly conduct the analysis requested.

To Submit Sample:

1. Once the required amount of groundwater has been collected into the proper sampling container, seal the container tightly with a screw cap lid.
2. Properly affix a label with the sample name, date and analysis.
3. Be sure to fill out the Chain of Custody (COC) form correctly and accurately and ship it along with the samples. A COC form is required for QA/QC purposes.
4. Once the bottles have been correctly labeled, place them in the designated cooler. Be sure to fill the remaining space in the cooler with blue ice or regular ice that has been double bagged in Ziploc bags to keep the entire shipment around 4°C, especially during the summer months.
5. All paperwork to be sent with the samples should be placed within a waterproof pouch or Ziploc bag and placed on top of the samples or affixed to the inside lid of the cooler.
6. Seal the cooler lid with a strong packaging tape.

SHIPPING INSTRUCTIONS

Packaging Samples:

1. Samples should be shipped in a cooler with ice or blue ice for next day delivery. If regular ice is used, the ice should be double bagged.
2. A chain of custody form must be included with each shipment of samples. Access our chain of custody at www.microbe.com.

Shipment for Weekday Delivery:

Samples for weekday delivery should be shipped to:

Attn: John Stephen
Microbial Insights (Australia) Pty Ltd
C/O AGRF Ltd
Plant Genomics Centre
Hartley Grove
Urrbrae
SA 5064

**Please note Saturday delivery is not available.